

ECAL DCS Beginners Manual

authors

June 3, 2020

Contents

1	Introduction	3
2	How to access the DCS	4
3	The ECAL DCS panels	4
3.1	Supervisor panel	5
3.2	ESS Panel	7
3.3	Low Voltage (LV) and High Voltage Panels	7
3.4	AirTemp, PTHM, Cooling Panels	8
3.5	Preshower Panels	9
4	Detector Partition	10

1 Introduction

The target of this manual are beginners who wants to approach for the first time the CMS ECAL DCS. The purpose of this manual is to explain its very basics. Throughout the whole document, images of the DCS will be showed and described.

The pictures used are taken from a virtual machine, that simulates the DCS. For this reason the status of some components is not the ideal one, as it will be explained in detail in every occasion.

2 How to access the DCS

(Explain in detail how to access the DCS) The access to the CMS ECAL DCS can be done remotely. It is important to have on your laptop a client of remote desktop connections. For example for Ubuntu there is the Remmina client. If you are not physically at CERN, you must do a preliminary step. You should access with the remote desktop client to the server "cernts.cern.ch", which divides the CERN Network to the outside. If you are at CERN you can skip this passage, because you don't need to access the CERN Network.

Once inside the CERN Network, both if you have used the client or if you already are at CERN, you should access to the server "ecal-dcs-vir06" (HERE we should insert the real server).

3 The ECAL DCS panels

Once entered in the DCS, the first screen that can be seen in the one in Fig. 1. The first thing to do is to log in with your personal account, with a right click on "NO USER", written in the upper right corner of the panel, as suggested by the red arrow in 1, or clicking on the close key.




Figure 1: Principal window of the DCS. The red arrow indicates where to click to login, as explained in the text.

Consequently the window showed in Fig. 2 on the left is going to appear. There you should insert your username and your personal password. If everything is correct, there will not be any error and the right corner of the DCS panel will appear as in Fig. 2 on the right, with your username. The last action to do to complete the login is to choose your role. You can right click where indicated by the red arrow in Fig. 2 on the right and choose your role, that most probably is "ECAL_Supervisor_Operator"




Figure 2: *On the left.* The log in window that pops up when clicking on the key on the main panel. You should insert your username and password to log in. *On the right.* After the log in, your username will show in the main panel. The red arrow indicated where to right click to change your role.

3.1 Supervisor panel

The principal panel in the one that can be seen as soon as you enter in the DCS. In Fig 3 three important areas of the main panel have been highlighted, and they will be explained in the following.

1. The first area includes different tabs. Each tab represents a different ECAL subsystem, which status must be monitored. The first, the "Supervisor", is the main panel, that is a summary of all the others. The details of each panel will be deeply explained in the following sections (3.2 - 3.5).
2. In this area it is possible to carry out a check of the ECAL status considering its different partitions. Each partition has a "tree" composition, because it is possible to access all the different subparts that compose the partitions. This important feature, that will be explained better in Section ??, is very useful because, if there is a problem, it is possible to find exactly which component has this problem.
3. The third area is the summary of the status of each partition of the ECAL. The graphical representation is very useful to have a fast and easily readable summary of the ECAL status.

In the upper part there are the Supermodules of the ECAL Barrel (EB), plus and minus. Below these there is the representation of the ECAL Endcap (EE), divided into Near (N) and Far (F), and of the Preshower (ES). The latter is attached to the Endcaps. The legend on the bottom helps in understanding the meaning of each color, which represents a different status of the calorimeter:

- ON → when the partition is working;




Figure 3: Main panel of the DCS. Three important areas of the panel have been circled and are explained in the text.

- OFF, STANDBY → when the partition is off or in standby;
- ERROR → the partition with this color has an error;
- PARTLY_ON, LV_ON_HV_OFF → the high voltage is OFF and the low voltage is ON;
- OFF_LOCKED → the ECAL is OFF in a safety mode;
- NO_DATA → the ECAL is not acquiring data.

From Fig. 3 it is possible to see that everything is in status ON, except for the ES, which is in state ERROR. As mentioned in the Introduction, all the screens of the DCS showed in this Manual are taken from a Virtual Machine. In this particular virtual machine the

ES components were not working, but usually, in the actual DCS, there should be no errors.

3.2 ESS Panel

Supervisor	ESS	LV	HV	Air Temp	PTHM	Cooling	ES LV	ES BV	ES Env	ES Heating		
Safety Conditions for						ECAL Barrel Plus	DI	DSS SHUTO.	POWER OK	COOL W DOG	MAGNET OFF	HEARTBEAT STABLE
Overheat.		Critical Ov.		Water leak		Cooling		LV C.Flow		LV C.Leak		
EBP01	Ok	Ok	<disabled>		OK		X3J01 (OK)		X3J01 (OK)			
EBP02	Ok	Ok	<disabled>		OK		X3J01 (OK)		X3J01 (OK)			
EBP03	Ok	Ok	<disabled>		OK		X3J01 (OK)		X3J01 (OK)			
EBP04	Ok	Ok	<disabled>		OK		X3J01 (OK)		X3J01 (OK)			
EBP05	Ok	Ok	<disabled>		OK		X3J01 (OK)		X3J01 (OK)			
EBP06	Ok	Ok	<disabled>		OK		X3V01 (OK)		X3V01 (OK)			
EBP07	Ok	Ok	<disabled>		OK		X3A01 (OK)		X3A01 (OK)			
EBP08	Ok	Ok	<disabled>		OK		X3A01 (OK)		X3A01 (OK)			
EBP09	Ok	Ok	<disabled>		OK		X3A01 (OK)		X3A01 (OK)			
EBP10	Ok	Ok	<disabled>		OK		X3A01 (OK)		X3A01 (OK)			
EBP11	Ok	Ok	<disabled>		OK		X3A01 (OK)		X3A01 (OK)			
EBP12	Ok	Ok	<disabled>		OK		X2R02 (OK)		X2R02 (OK)			
EBP13	Ok	Ok	<disabled>		OK		X0R11 (OK)		X0R11 (OK)			
EBP14	Ok	Ok	<disabled>		OK		X0R11 (OK)		X0R11 (OK)			
EBP15	Ok	Ok	<disabled>		OK		X0R11 (OK)		X0R11 (OK)			
EBP16	Ok	Ok	<disabled>		OK		X0U11 (OK)		X0U11 (OK)			
EBP17	Ok	Ok	<disabled>		OK		X2U01 (OK)		X2U01 (OK)			
EBP18	Ok	Ok	<disabled>		OK		X2U01 (OK)		X2U01 (OK)			
Description of safety condition notifications												
OVERHEATING	SAFETY PROBLEM is currently present!					OVERHEATING	NO safety problem. Interlocks are still set.					
OK	NO safety problem. Interlocks are not set.					<disabled>	Condition is EXCLUDED from the ESS action.					

Figure 4: ESS Panel

3.3 Low Voltage (LV) and High Voltage Panels

These two panels can be accessed clicking in the tabs LV and HV respectively, showed in Fig. 3 in the first circled area. As can be checked in Fig 5, they are very similar to the principal panel, the *Supervisor* tab. The difference is that the principal panel shows a summary of the status of the ECAI and the Preshower, while here there is only information about LV (and

HV) in the ECAL. Indeed there is no information about the ES, which can be found in other tabs (see section 3.5).




Figure 5: *On the left.* Low Voltage panel. *On the right.* High Voltage panel. They show the status of the LV (HV) in the ECAL Barrel and in the ECAL Endcap, using the same color legend as in the Supervisor panel.

Let's now consider the case in which the HV in the positive Barrel is OFF. In this situation is important to understand the role of each tab to really understand what is happening. In the Supervisor panel (Fig. 6), where we have a summary of all the ECAL situation, there is the EB+ in the yellow status, which means LV_ON_HV_OFF.

The status can be also checked in the LV and HV panels, as shown in Fig. 7. On the left the LV panel is shown, and the status is all ON, while the figure on the right, that represents the HV panel, shows that the EB+ is OFF.

In conclusion, depending on the tab that we are visiting, it is possible to access to a detailed or summary information.

3.4 AirTemp, PTHM, Cooling Panels

The parameters of Air Temperature (AirTemp), of Pressure Temperature Humidity Monitoring (PTHM) and Cooling are shown in the respective tabs. This information comes from the Safety System ECAL. These are only control parameters, they can not be modified, but they need to be checked to be sure that everything is working. The screens in Fig. 11 and Fig. 9 show how these panels appear. Again they are very similar to the principal panel. The PTHM and Cooling panels also show the temperature of the water, in place of the violet boxes, which represent the NO_DATA state. The figures here reported are taken from a Virtual Machine, so the Safety System outputs are not connected, as it is possible to be seen in Fig. 9.

3.5 Preshower Panels

All the other panels are dedicated to the Preshower.




Figure 6: This is the Supervisor panel of the DCS with the EB+ with the HV OFF. This can be seen by the graphical representation of the ECA, in which the EB+ is yellow, which means PARTLY_ON, and in the left side of the window, where it is explicitly written, in correspondence of "Barrel Plus", that it is in the status of LV_ON_HV_OFF.




Figure 7: These are the screens of the ECAL DCS panels LV and HV in the situation described in Fig. 7. *On the left.* The LV panel, in which everything is ON. *On the right.* The HV panel, in which the EB+ barrel has the OFF status.

4 Detector Partition

The second area showed in Fig. 3. (Add an error to one component and show that you can find the exact part with the error using the tree composition of this part of the main panel)




Figure 8: *On the left.* The AirTemp panel. *On the right.* The PTMH panel. The peculiarity here is the information about the water temperature, in place of the violet boxes, that here is not showed because the virtual machine from which the pictures are taken is not connected to any Safety System.




Figure 9: The Cooling Panel.




Figure 10: *On the left.* The AirTemp panel. *On the right.* The PTMH panel. The peculiarity here is the information about the water temperature, in place of the violet boxes, that here is not showed because the virtual machine from which the pictures are taken is not connected to any Safety System.




Figure 11: *On the left.* The AirTemp panel. *On the right.* The PTMH panel. The peculiarity here is the information about the water temperature, in place of the violet boxes, that here is not showed because the virtual machine from which the pictures are taken is not connected to any Safety System.